SOLLUME ESTHÉ

LED REAL MASK

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Product name	SOLLUME ESTHÉ LED REAL MASK		
Warranty period	1 year from the purchase date		
Purchase date	(year)	(month)	(date)
	* Please write down the purchase date. Warranty might be invalid when the date is not written.		
Customer name		Contact number	
Purchase location			

Free warranty is provided under warranty regulations if failure occur when the product is used normally under precautions such as instructions within the warranty period. If warranty regulations apply to service at cost, you must pay for the service. For warranty service, contact [G&J] customer service center with the warranty after checking warranty regulations.

※ Customer service center : 070-8850-0111

[Warranty regulations]

The object of warranty is the product body.

[Free service]

Free warranty is provided if failure occur when the product is used normally within 1 year from the purchase date (warranty period).

When the product is used for business purposes, free service is not provided.

[Service at cost]

Service for the following cases is provided at a cost even under the warranty period:

1. When failure occurs due to negligence of a consumer (fall, water damage, breakage, damage after purchase);
2. External changes occurred during use (scratch, deformation, discoloration, etc.);
3. When the product is disassembled and modified by a consumer arbitrarily;
4. When damage is caused by using components other than parts included in this product or not following 'Cautions'; and
5. When the warranty is not provided, the purchase date is not written on the warranty, or contents are modified intentionally.

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▶ Product care

- 1) Wipe the part contacted with skin using a cotton pad, soaked in toner, or a dry towel (tissue, etc.) after using the product.
- 2) Do not wash with water or other liquids as the product is not waterproof.
- 3) Do not wash with chemical agents such as thinner and benzene.

▶ Storage

- 1) Store the mask straight up.
- ※ Caution : Mask and remote controller are not separated. Do not pull the connecting cable hard.
- 2) Store the product at a room temperature.
- 3) Avoid direct sunlight and do not store in a humid place.
- 4) Do not store with sharp objects that may damage the product.

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1. Read instructions thoroughly and understand how to use the product before use.
2. Do not use the product for any other purpose than intended use.
3. Do not open the frame from side to side with too much pressure.
4. Do not store in a humid place or drop the product in water.
5. Do not use on your face if you have wounds or dermatitis, and follow directions of a doctor
6. If abnormal reactions such as itchiness or rash occur during use, stop using the product immediately and follow directions of a doctor.
7. Do not wash with water as the product is not waterproof.
8. Do not clean the product with chemical agents or detergents.
9. Do not beam LED directly on eyes.
10. Do not use the product thermal products (the product may be deformed or discolored).
11. Charge the battery frequently (within 1 month) when the product is unused for a long time.
12. Exchange or refund is not available when failure is occurred due to negligence of a consumer.
13. Free service is provided within 1 year from purchase (warranty period) (consumer negligence is excluded).

※ Consumer negligence

- Failure caused by negligence of a consumer (fall, water damage, shock, breakage, forced pressure, etc.)
- Failure caused by repair of a person other than appointed service technician or partner technician.
- Failure caused by indentation or mistake of a consumer.
- Failure or damage caused by using components or parts that are not approved.